



Endurance™

Tilt Chair

*Unparalleled strength & style.
Complement your Endurance™ Instrument Stand.*

Endurance™ Tilt Chair



The Reichert® **Endurance™ Tilt Chair** meets the needs of the eye care professional who demands the very best in comfort, style, functionality and durability.

Endurance™ Tilt Chair Features:

- **Premium quality.** Heavy-duty construction provides long-lasting performance.
- **Effortless maneuverability.** Counterbalanced 45° tilt and 360° rotation enables any staff member to easily position the patient as needed for different procedures and testing. Raise and lower the chair using the foot pedal, controls on the slit lamp arm or the stand control panel.
- **No more fluid leaks.** The industrial strength, Swiss designed linear actuator lift mechanism is quiet, durable and maintenance free, with a lift capacity of 500 lbs.
- **Easy access.** The stainless steel footrest and upholstered arms rotate upward and out of the way for convenient patient ingress and egress.
- **Maximum patient comfort.** Deep cushions offer a generous seat width, while premium upholstery remains soft to the touch without sacrificing durability.
- **Color Options.** Color options are available to complement your practice. See Endurance Tilt Chair at www.reichert.com/endurance for more details.
- **Peace of mind.** 2 year standard warranty on the chair, 5 year warranty on linear actuator lift mechanism.



Reichert® **Endurance™**
Tilt Chair & Stand

Endurance™ Tilt Chair Specifications

Catalog Number	15250	Chair Rotation	360°
Base Dimensions	20 in. W x 22 in. D (52 x 56 cm)	Chair Tilt	45°
Seat Dimensions	21 in. W x 19 in. D (53 x 48 cm)	Weight	270 lbs. (123 kg)
Depth (Fully Tilted)	50.25 in. (128 cm)	Fuses	6 Amps Main
Height to Headrest	50.5 in. (128 cm)	Electrical	120 VAC, 60 Hz
Chair Elevation	22.5 to 30.5 in. (57 to 78 cm)		



Visit www.reichert.com or contact your Authorized Reichert Technologies Distributor.

Important Notice: For your own protection, purchase instruments from and have them serviced only by Reichert Technologies or an Authorized Reichert Technologies Distributor.

15250-110-Rev-C ©2013 AMETEK, Inc. All rights reserved. - Specifications subject to change without notice. Reichert and Reichert Technologies are registered trademarks of Reichert, Inc. Endurance and "Advancing Eye Care. Preserving Sight." are trademarks of Reichert, Inc. AMETEK is a registered trademark of AMETEK, Inc.



Reichert, Inc.
3362 Walden Avenue, Depew, NY 14043
Toll Free: 888-849-8955 • Tel: (716) 686-4500
Fax: (716) 686-4545 • www.reichert.com

ISO 9001 Certified
Reichert products are designed and manufactured under quality processes meeting ISO 9001 requirements.