

500

Series Hydraulic Stools



800.735.0357
haag-streit-usa.com



Model 540 Hydraulic Stool with 15" Round Seat



Model 546 with Tilting Adjustable Backrest



Model 548 with Waterfall Seat



Model 556 with Wrap-around Arm

500

Series Surgical Stools



800.735.0357
haag-streit-usa.com

500 Series Surgical Stools

The base exceeds the circumference of the stool to provide stability even in the sit-stand position.

Swivel of stool seat is locked by lifting floor pedal. A convenient foot pedal easily raises or lowers the stool height without using hands.

Choice of Three Heights

Standard: 17½" to 25"

Medium: 20½" to 28"

High: 24" to 31½"

Finished for Durability

- Select from 14, easy-to-maintain vinyl colors. Other vinyl colors available by request.
- All exposed metal with exception of the aluminum base is chrome-plated steel.
- Base and legs finished in tough, black epoxy enamel.

Optional Foot Ring

Order a foot ring with your stool, or at a later date.

Exceptional Stability

All 500 Series stools have five legs to evenly distribute weight and ensure stability.

Tech Specs

Base diameter 19¾"; footrest: 16" diameter; 2" spacer available to elevate ring 2" to 4".

Shipping weights: Model 540 (47 lbs.), Model 546 (51 lbs.), Model 548 (57 lbs.), Model 556 (63 lbs.)

	540	546	548	556
Seat Adjustment	Foot Pedal	Foot Pedal	Foot Pedal	Foot Pedal
Seat	15" circular	15" circular	14 x 16" rectangular	15" circular
Backrest – tilt, horizontal and vertical adjustment	None	Yes	Yes	Yes
Wrap-around arm	No	No	No	Yes
Seat	2" polyfoam	2" polyfoam	2" polyfoam	2" polyfoam
Five legs for extra stability	Yes	Yes	Yes	Yes
Foot pedal	Yes	Yes	Yes	Yes
Floor lock	Optional	Optional	Optional	Optional