

BABY SCALE

PRODUCT NUMBER 8310





BABY SCALE

PRODUCT NUMBER 8310

Features

- > Beautifully designed handy baby scale made of easy to clean plastic
- > Removable baby weighing tray: can be used as personal scale for weighing children up to 20kg
- > Quattro-Tronic: fully electronic scale with 4 load cells for rapid precise weight measurement
- > No wear parts and maintenance-free
- > Easy to read LCD with 20mm high digits
- > Simple one-key operation for on/off and tare

- > Weight lock function: weight display remains stable if baby moves about
- > Blankets can be tared
- > Lightweight sturdy carry case to protects the scale when transported
- Battery power supply with 9V block battery and energy saving automatic power-off

Ordering information

Order Number	Order Number	Weighing Range	Divisions	Dimensions	Intrinsic weight
	Carry case	kg	g	mm	kg
8310.00.090	5030.00.071	20	10	530 X 240 X 100	4

Your Soehnle Professional dealer:

www.soehnle-professional.com