


Are high quality slit lamp images important to you?

EyeCap™ SL

A COMPREHENSIVE, EASY TO USE
IMAGE CAPTURE AND MANAGEMENT
SYSTEM FOR THE CAPTURE OF HIGH
QUALITY SLIT LAMP IMAGES.



Contact your
HAAG-STREIT dealer or



EyeCap™ SL


- High quality images
- Simple operation
- Live on screen view
- 3CCD, 1 CCD & digital capture
- Vast on-line storage
- Movie function*
- DICOM export

EyeCap™ SL is a complete image capture and management application that has been developed with cooperation from several leading ophthalmic clinicians, ophthalmic photographers and optometrists as a means to simply, accurately and securely document the eye via the slit lamp.

Images may be captured directly from slit lamps equipped with digital or video cameras as well as any other device with video or TWAIN output. The live image can be viewed on the computer screen and is then instantly stored in the database at the press of a footswitch. Segments of real time video may also be captured providing practitioners with an excellent tool for recording and monitoring of ocular pathology.

EyeCap™ employs a unique image archiving systems that keeps a copy of all images captured on-line providing the ability to store several hundred thousand images using conventional computer hardware. This intelligent method of image storage provides clinicians with fast-track access to images from any patient visit regardless of when the original images were acquired. All images are stored in a secure environment but can be easily exported for referral or teaching and can be easily used with other applications and third party software.

* With Sony IEE 1394 (Firewire) Camera, check with your distributor for current recommended model.



Images courtesy of Leyenburg Hospital, The Netherlands